

SPECIALTY
APPLICATIONS

CT Dental Viewer

Neurology

Dentistry

PANORAMIC VIEWS OF THE JAW





With a single CT volume that captures the upper or lower jaw, define a pathway to create panoramic reconstructions that include automatically generated cross-cut images and volume-rendered views. Make measurements that are represented on 3D volume-rendered views.

Cardiovascular

Respiratory

Key features:

- Create panoramic images simply by tracing the dental arch
- View cross-cut from any direction
- Visualize and measure the relationship between the mandibular canal and adjacent teeth
- Create a cine output of serial cross-cut images



Japan Office

1-4-28 Mita, Minato-ku, Tokyo, 1080073, Japan
Tel: +81 3 5427 1903 <http://www.zio.co.jp>

US Office

1301 Shoreway Road, #325, Belmont, CA 94002-4105, USA
Tel: +1 650 413 1300 <http://www.ziosoftinc.com>