




PET/CT Fusion

FAST FUSION OF FUNCTIONAL AND ANATOMIC DATA

The PET/CT Fusion protocol provides flexible and quick image overlay analysis. Standard Uptake Value (SUV) calculations are enabled with regions and volumes of interest (ROIs and VOIs).

Key features:

- Registration controls for any position and size differences between datasets
- View multiple images and slice axes simultaneously
- Multiple ROIs yield SUV data from different regions simultaneously
- Open PET, MRI or NM functional data with CT or MRI anatomical data for fusion



Japan Office

1-4-28 Mita, Minato-ku, Tokyo, 1080073, Japan
 Tel: +81 3 5427 1903 <http://www.zio.co.jp>

US Office

1301 Shoreway Road, #325, Belmont, CA 94002-4105, USA
 Tel: +1 650 413 1300 <http://www.ziosoftinc.com>